ProteoSil

Instruction Manual

1 Introduction

Thank you for purchasing this product.
Please read this instruction manual carefully and use the product properly to ensure the best performance.

2 Handling Precautions

2-1 Precautions

- Do not drop or hit the column. Strong impacts may cause column degradation.
- The column provides high pressure resistance because it is packed using a high-pressure slurry method. However, in order to ensure long-term stable use, we recommend using the column at pressures not exceeding than those shown in the table below.

Particle diameter	Column size	Recommended pressure
1.9 μm	All sizes	80 MPa or less
3 μm HP	All sizes	50 MPa or less
3 - 10 μm	I.D. 0.3 to 50 mm	20 MPa or less

- Be careful of sudden pressure fluctuations.
 - When removing the column, wait until the pressure gauge indicates 0.
 - Please note that slow operation of the sample injection valve will cause sudden pressure changes at the column inlet.
- To the extent possible, dissolve the sample in a solvent of the same composition as the eluent (initial solvent for gradients).
Injecting a large amount of sample dissolved in a solvent with a higher solubility than the eluent can result in reduced separating capacity and precipitation of the sample at the inlet of the column.
- To replace eluents that do not mix with each other (e.g., changing from a normal-phase type to reversed-phase type), an intermediate-polar solvent such as 2-propanol should be poured in between.
At such a time, since alcohol-based solvents are subject to high pressures during feeding, adjust the flow rate while paying attention to the upper pressure limit of the column, and feed the solvent at a rate of at least 10 times the column capacity.
- If a peak that elutes quickly tails off, the cause may be due to dead volume. Check if the connection piping at the column joint is inserted all the way to the end.
For piping to the injector and detector, select piping with an inner diameter and length suitable for the applicable column inner diameter and corresponding analytical system. The effect of piping is particularly significant when analyzing at low flow rates using semi-micro columns, etc.
- Pressure rise and peak cracking may be caused by clogging or contamination of the filter at the column inlet.
 - Filter the eluent through a membrane filter of 0.45 μm or less before use.
 - Filter the sample solution with a syringe filter of 0.45 μm or less before injecting.
 - Column clogging can be prevented by using Cartridge E or Cartridge E 2.1x10 mm UHPLC.
- Ensure that the column is well balanced with the eluent before use.
(Reversed-phase mode using HILIC or ion-pair reagents may take longer to stabilize.)
- Depending on LC/MS analysis conditions and the type of equipment, precipitates may be seen at the column outlet. Check the effect of precipitation before use.


- Refer to the table below for the applicable pH range and upper temperature limit of the column. (Refer to the GL Sciences website for the applicable pH range and upper temperature limit for columns not listed in the table.)
- Please note that exceeding the upper temperature limit can lead to premature deterioration.

Product name	Applicable pH range (Normal 20 to 40 °C)	Upper limit temperature	
		60°C (pH 1-7)	50°C (pH 1-10)
ProteoSil 200-C18, ProteoSil 100-C18	1-10 *1,2,3	60°C (pH 1-7)	50°C (pH 1-10)
ProteoSil 200-C8	1-10 *1,2,3	60°C (pH 2-7)	50°C (pH 1-9)
ProteoSil HILIC	2-8.5 *2	60°C (pH 2-7)	50°C (pH 2-8.5)
ProteoSil 300-C18 / 300-C8 / 300-C4 / 100-SEC / 300-SEC	2-7.5 *2	60°C (pH 2-7)	50°C (pH 2-7.5)

- *1 Column life varies greatly depending on conditions such as applicable pH, temperature, and eluent composition. To ensure stable and long-term use, we recommend lowering the column temperature, using low concentrations of buffer salts and additives, and analyzing with eluents containing organic solvents. For analysis at pH 1-2, we recommend using TFA, formic acid, acetic acid, phosphate, etc.
For analysis at pH 9-10 within the applicable pH range, we recommend the use of an organic buffer solution of about 5 mM (triethylamine, etc.). Use at pH 2-8 when analyzing only buffer solutions that do not contain organic solvents.
- *2 In order to prevent premature deterioration, do not to exceed the pH of the eluent described above.
- *3 For use at pH 1-2 or pH 9-10 within the applicable pH range, we recommend analysis at low temperatures and use of eluents containing organic solvents such as methanol.

2-2 UHPLC PEEK / PEEK Column Handling Precautions

- All wetted parts of the UHPLC PEEK / PEEK column are made of resin. The column port in particular has a very delicate structure. Therefore, if the column is handled in the same manner as a stainless-steel column, the wetted parts may become damaged.
- Since it can be connected with less torque than a stainless steel column, when connecting to equipment piping, etc., make sure that there is no leakage of fluid, and tighten with sufficient torque (about 0.8 N-m) to prevent leakage.
- Do not use a male nut with a deformed tip when using it for prolonged periods of time or repeatedly. Since the column port may become damaged, we recommend using new resin male nuts when making connections.
[Recommended fitting: 2-piece type resin male nut.]
- Feeding THF or chloroform through it for a prolonged period of time can lead to premature column degradation.
- When installing or removing the column, use the connecting male nut and column end nut portion.

3 Product Features

This product is manufactured under our strict standards at each production process for matrix silica gel, chemical modification, and column performance inspection after packing. This ensures that the product can be used at a consistent level of quality with peace of mind.

4 Confirmation of Contents

- Check for any abnormalities in column appearance, packaging, etc.
- Check for any errors in packing material name, column size, etc.
- Check to see that the column performance report is included.
The column performance report includes information such as the packing material lot No., column serial No., and column performance test eluent. Please keep it in a safe location.
- In principle, columns are sealed with the eluents described in the column performance report.

5 Specifications

Type	Main column	Eluent sealed at shipping
Normal-phase type	ProteoSil 100-SEC / 300-SEC	Hexane / ethanol Mixed eluent
Reversed-phase type	ProteoSil 100-C18 / 200-C18 / 300-C18 / 300-C8 / 300-C4 / 200-C8	Acetonitrile / water Mixed eluent
HILIC type	ProteoSil HILIC	Acetonitrile / water Mixed eluent

6 Storage

- If an eluent containing a buffer salt or ion-pair reagent is used in a reversed-phase column, wash it thoroughly with an eluent free of salts.
- Store the reversed-phase column after replacement with an organic solvent such as acetonitrile or methanol.
- When using a reversed-phase eluent, wash ProteoSil 100-SEC or 300-SEC columns with water/acetonitrile in a 50/50 ratio, etc. Store after replacement with 100% acetonitrile.
- Wash the normal-phase column with ethanol or 2-propanol.
Since alcohol-based solvents are, generally speaking, subject to high pressures during feeding, pay attention to the upper pressure limit of the column and feed it at a reduced flow rate if necessary. Store after replacement with 100% hexane.
- After using a HILIC type column with an eluent containing buffer salts, feed it with an eluent with a high water concentration (50% water) to remove hydrophilic substances, and store it after replacement with a high acetonitrile concentration (80% or higher acetonitrile).
- When storing the column, tightly plug it with the provided plug and store it in a cool place where there is little temperature change or humidity.
*For UHPLC PEEK / PEEK columns, mount the plug in the same manner as when connecting to equipment, while being careful not to overtighten it.

ProteoSil Series products are manufactured, inspected, packaged, and shipped under strict quality control, but in the unlikely event of a defect, please contact one of our branches, sales offices, or distributors.

However, we do not provide any warranty against damage, problems related to service life, or deterioration caused by use of the product in a manner not in compliance with this instruction manual.

The specification and appearance are subject to change without notice.

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