REFERENCE GUIDE

# Optima series ICP-OES

# Sample Introduction Components



#### **NEBULIZERS**

Suitable for Optima 2X00, 4X00, 5X00, 7X00 DV and V

Description	Part No.
SeaSpray <sup>™</sup> U-Series <sup>™</sup> Nebulizer, 2 mL/min uptake	N0775345
For best performance or for samples with high salt content	
Conikal™ U-Series™ Nebulizer, 2 mL/min uptake	N0777487
For routine analysis	
MicroMist™ U-Series™ Nebulizer, 0.6 mL/min uptake	N0775346
MicroMist™ U-Series™ Nebulizer, 0.2 mL/min uptake	N0775347
MicroMist™ U-Series™ Nebulizer, 0.1 mL/min uptake	N0775348
MicroMist™ U-Series™ Nebulizer, 0.05 mL/min uptake	N0775349
For organics or for samples with limited volume	
Slurry™ U-Series™ Nebulizer, 4 mL/min uptake	N0777488
For slurries or suspensions – ideal for wear metals in oil	mesonapper to
PolyCon™ Nebulizer, 2 mL/min uptake	N0777490
For samples containing HF	
VeeSpray <sup>™</sup> Ceramic Nebulizer	N0777489
Most versatile nebulizer – handles samples containing particulates and high salt content and HF	



# **SPRAY CHAMBERS**

# Suitable for Optima 2X00, 4X00, 5X00, 7X00 DV

Description	Part No.
Tracey™ Cyclonic Spray Chamber	N0776052
Tracey™ Cyclonic Spray Chamber with Helix™	N0775351
For best sensitivity	
Twister™ Cyclonic Spray Chamber (7 mm baffle)	N0776053
Twister™ Cyclonic Spray Chamber with Helix™ (7 mm baffle)	N0775352
For best precision	
Twister™ Cyclonic Spray Chamber (4 mm baffle)	N0776090
Twister™ Cyclonic Spray Chamber with Helix™ (4 mm baffle)	N0775355
For organics	
Jacketed Twister™ Cyclonic Spray Chamber with Helix™	N0774101
For organics – chiller required	
Cinnabar™ Cyclonic Spray Chamber with Helix™	N0777494
Twinnabar™ Baffled Cyclonic Spray Chamber with Helix™	N0777495
For samples with limited volume	
Tracey™ TFE Cyclonic Spray Chamber with Helix™	N0777496
For samples containing HF	

# **SPRAY CHAMBERS**

# Suitable for Optima 4X00, 5X00, 7X00 V

Description		Part No.
Twister™ Baffled Cyclonic Spray Chamber		N0771501
Twister™ Baffled Cyclonic Spray Chamber with Helix™		N0775356
For aqueous or organics		
Tracey™ TFE Cyclonic Spray Chamber with Helix™	The same of the sa	N0777497
TC-44 Cyclonic Spray Chamber		N0771550
For samples containing HF		

# **TORCHES AND TORCH PARTS**

# Suitable for Optima 2X00, 4X00, 5X00, 7X00 DV

Description	Part No.
D-Torch™ for Optima™ 2/4/5/7000 DV	N0777508
Demountable torch reduces torch replacement costs	
Base and Inner Tube for D-Torch™ Optima 2/4/5/7000 DV	N0777509
Quartz Outer Tube for D-Torch™ Optima 2/4/5/7000 DV	N0777510
Safety Gloves – Medium	N0777421

# **TORCHES AND TORCH PARTS**

Suitable for Optima 4X00, 5X00, 7X00 V

Description	Part No.
D-Torch™ for Optima 4/5/7000 V	N0777051
Demountable torch reduces torch replacement costs	
Base and Inner Tube for Optima 4/5/7000 V D-Torch™	N0777054
Quartz Outer Tube for Optima 4/5/7000 V D-Torch™	N0777053
Ceramic Outer Tube for Optima 4/5/7000 V D-Torch™	N0777052
For extended life – ideal for oil analysis	
Retaining Ring for D-Torch™	N0777055
Safety Gloves – Medium	N0777421
Torch Adaptor with Auxiliary Port	N0774070

# **INTERNAL STANDARD KIT AND SPARES**

# Suitable for all models

Description	Part No.
Internal Standard Kit with orange/blue flared end tubing	N0774068
HF-resistant Internal Standard Kit with orange/blue flared end tubing	N0774067
EzyFit™ glass mixing T-piece	N0777461
HF-resistant mixing T-piece	N0777466
EzyFit™ sample tube, twin (pack of 6)	N0777462
UniFit™-EzyFit™ sample tube (pack of 6)	N0777061
EzyFit™ sample tube, 0.5 mm ID (pack of 10)	N0777463
EzyFit™ sample tube, 0.25 mm ID (pack of 10)	N0777464
Internal standard probe	N0777465
Flared end orange/blue 2 tag PVC peristaltic pump tubing 0.25 mm ID (pack of 6)	N0773112

# **CONNECTORS AND MISC. PARTS**

Description	Part No.
Hose clip (pack of 10)	N0773197
UniFit™ Connector with 0.75mm ID x 1.3 mm OD x 700 mm long sample tube (pack of 10)	N0777512
UniFit™ Connector with 0.50 mm ID x 1.3 mm OD x 700 mm long sample tube (pack of 10)	N0777513
UniFit™ Connector with 0.25 mm ID x 1.3 mm OD x 700 mm long sample tube (pack of 10)	N0777514
EzyFit™ Connector with 0.75 mm ID x 1.3 mm OD x 700 mm long sample tube (pack of 10)	N0777435
EzyFit™ Connector with 0.50 mm ID x 1.3 mm OD x 700 mm long sample tube (pack of 10)	N0777436
EzyFit™ Connector with 0.25 mm ID x 1.3 mm OD x 700 mm long sample tube (pack of 10)	N0777437
UniFit <sup>™</sup> Drain Connector with 1.3 mm OD X 0.75 mm ID x 700 mm long drain tube (pack of 3)	N0774077
UniFit™ Drain Connector with 1/16 inch OD X 0.75 mm ID x 700 mm long drain tube (pack of 3)	N0774080
EzyLok™ connector for 1/4 inch tubing (pack of 4)	N0774100
EzyLok™ Kit	N0777413
Nebulizer Adaptor for Cyclonic Spray Chamber 16/6	N0776006
O'ring Set for Nebulizer Adaptor 16/6	N0776007
Helix™ Locking Screw	N0777438
Helix™ Seal (pack of 4)	N0777439

#### **ACCESSORIES**

Description Part No.

Argon Humidifier

For samples with high salt content

N0775331

Eluo™ Nebulizer Cleaner To unblock nebulizers



N0775358

IsoMist™ Kit for Optima 4300/5300/7300 V IsoMist™ Kit for Optima 2/4/5/7000 DV Programmable temperature spray chamber Improves stability, simplifies analysis of volatile organics



PerkinElmer, Inc. 940 Winter Street Waltham, MA 02451 USA P: (800) 762-4000 or (+1) 203-925-4602 www.perkinelmer.com

