

Agilent Triple Quadrupole LC/MS Supplies

Make sure that you don't run out of essential columns and supplies

Use this handy list to keep your shelves stocked—and keep your lab productive.



Agilent Ultivo LC/TQ	Part Number
Ion injector capillary	
Ion injector	G6301-80004
Narrow-bore capillary cap	G1960-20310
Desolvation assembly, ion optics, and vacuum manifold	
Canted coil spring (4)	G2571-67001
High-gain electron multiplier horn	G2571-80103
Skimmer, 2 mm	G1969-20302
AJS spray shield, small	G1958-20050
Endcap assembly, multimode	G1978-60060
Calibrant delivery system	
50 mL CDS bottle + cap	G6301-80019
Rough pump (MS40+)	
Nitrogen filter for main gas supply (2)	RMSN-4
AFV 60 gold oil	6040-1444
Oil mist filter, MS40+	G1960-80039
Gas purifiers	
Nitrogen filter for main gas supply (2)	RMSN-4
Nitrogen gas filter for collision cell gas supply	RMSN-2

Agilent 6495A/B/C LC/TQ	Part Number
Ion injector capillary	
Multi bore capillary cap [G6495A/B]	G1964-20310
Capillary, 6 x 0.61 mm x 90 mm, 1G	G1964-80659
Capillary cap, 0.8mm ID [G6495C]	G1964-20028
Desolvation assembly, ion optics, and vacuum manifold	
Canted coil spring (4)	G2571-67001
High-gain electron multiplier horn	G2571-80103
AJS spray shield, multi bore	G1964-20050
Endcap assembly, multimode	G1978-60060
Calibrant delivery system	
Glass CDS bottle	9300-2576
Glass CDS bottle cap	9300-2575
Rough pump (MS120) [G6495C only]	
AVF 60 gold oil (2)	6040-1444
MS120 exhaust filter kit	X3702-68201
Rough pump (MS40+) [G6495A/B only, (2)]	
AVF 60 gold oil (2)	6040-1444
Oil mist filter, MS40+ (2)	G1960-80039

Agilent 6495A/B/C LC/TQ	Part Number
Gas purifiers	
Nitrogen filter for main gas supply (2)	RMSN-4
Nitrogen gas filter for collision cell gas supply	RMSN-2
Agilent G6407A/B and G6475A LC/TQ	
Ion injector capillary	
0.6 mm id fast switch capillary	G1960-80060
Narrow-bore capillary cap	G1960-20310
Ion Injector, 180 mm, 0.6 mm id, Dielectric	G7604-60000
Desolvation assembly, ion optics, and vacuum manifold	
Canted coil spring (4)	G2571-67001
High-gain electron multiplier horn	G2571-80103
Skimmer, 2 mm	G1969-20302
AJS spray shield, small	G1958-20050
Endcap assembly, multimode	G1978-60060
Calibrant delivery system	
Glass CDS bottle	9300-2576
Glass CDS bottle cap	9300-2575
Rough pump (MS40+)	
AFV 60 gold oil (1)	6040-1444
Oil mist filter, MS40+	G1960-80039
Gas purifiers	
Nitrogen filter for main gas supply (2)	RMSN-4
Nitrogen gas filter for collision cell gas supply	RMSN-2
Common-Use Parts	
Tools	
1/4 inch × 5/16 inch open-end wrench	8710-0510
1/2 inch × 9/16 inch open-end wrench	8710-0877
9/64 inch hex key	8710-2394
LC tools, long hex key kit (1.5, 2, 2.5, 3, 4 & 5mm)	5023-2524
10 mm hex key	8710-2612
4000 grit abrasive paper	8660-0827
Clean gloves	8650-0030
Clean chemical-resistant gloves	9300-1751
Clean lint-free cloth 15/pk	05980-60051

Common-Use Parts	Part Number
Cotton swab 100/pk	5080-5400
Long-nose pliers	8710-0004
Nebulizer adjustment kit	G1960-67470
T10 torx screwdriver	8710-1623
T20 torx screwdriver	8710-1615
Alconox powdered precision cleaner	5190-1401
Reagents	
LC/MS-grade or de-ionized (18 MΩ/cm) water (4 L)	8500-2236
HPLC flushing solvent	G1969-85026
LC/MS-grade acetonitrile (1 L)	G2453-85050
LC/MS-grade methanol (1 L)	5190-6896
Formic acid (5 mL)	G2453-85060
InfinityLab Ultrapure LC/MS solvent 1 L: Acetonitrile	5191-4496
InfinityLab Ultrapure LC/MS solvent 1 L: Methanol	5191-4497
InfinityLab Ultrapure LC/MS solvent 1 L: Water	5191-4498
Calibrant delivery system	
APCI-L tuning mix	G1969-85010
ESI-L tuning mix	G1969-85000
MMI-L tuning mix	G1969-85020
Rotor seal	0101-1409
CDS PEEK tubing	0100-1919
Ionization source	
ES nebulizer needle kit for AJS	G1958-60137
APCI nebulizer needle kit	G2428A
Corona needle	G1947-20029
Nebulizer assembly (ES)	G1958-65098
Nebulizer assembly (APCI)	G1946-67037
PM kits	
LC/MS PM kit, JetStream, 1 MS40+ RP	5190-9008
LC/MS PM kit, JetStream, 2 MS40+ RPs	5190-9009
LC/MS PM kit, MS120 RP	G2571-68003
Ultivo PM kit	G6301-68010

Agilent CrossLab services. From insight to outcome.

CrossLab is an Agilent capability that integrates services and consumables to support workflow success and important outcomes like improved productivity and operational efficiency. Through CrossLab, Agilent strives to provide insight in every interaction to help you achieve your goals.

CrossLab offers method optimization, flexible service plans, and training for all skill levels. We have many other products and services to help you manage your instruments and your lab for best performance.

Learn more about Agilent CrossLab, and see examples of insight that leads to great outcomes, at www.agilent.com/crosslab



Agilent
CrossLab

From Insight to Outcome

To place your order now, go to:
www.agilent.com/chem/lc-ms-supplies

DE.5424652778

This information is subject to change without notice.

© Agilent Technologies, Inc. 2023
Published in the USA, May 22, 2023
5990-8583EN

 **Agilent**
Trusted Answers